ENCODE DCC Antibody Validation Document


Date of Submission
Name: Email:
Lab
Lab
Antibody Name: Target:
Company/
Source:
Catalog Number, database ID, laboratory Lot Number
Antibody
Description:
Target
Target Description:
Species Toyant Species Hest
Species Target Species Host
Validation Method #1 Validation Method #2
Purification Polyclonal/
Method Monoclonal
Vicilia IIDI
Vendor URL:
eference (PI/
ublication
formation)
ease complete the following for antibodies to histone modifications:
your specifications are not listed in the drop-down box,
ease write-in the appropriate information
istone Name AA modified AA Position Modification
istone ivalile AA modified AA Fosition Iviodification

Validation #1 Analysis				
Insert Validation Image (click here)				

Antibody: TATA binding protein TBP Source: Abcam ab62126

Epitope: Mouse monoclonal antibody to amino acids 1-100 of human TBP

Validation 1: Immunoprecipitation (IP) in both CH12 and MEL cell lines



Arrow indicates immunoprecipitated band of expected size of TBP in both CH12 and MEL cell lines (~38 kDa). We used DMP crosslinker to covalently couple the antibody to the protein A beads in the IP experiment to get rid of antibody bands.

NE: nuclear lysate

S: supernatant after IP

IP: IP with tested antibody

IgG: IP with control IgG

Note: TBP might have phosphorylation or some other modification which caused the band shift

Validation #2 Analysis					
Insert Validation I	mage (Click here)				
Insert Validation Image (Click here)					